

SYMANSIS™
CELL SIGNALING SCIENCE

Easy50™ PI3K Inhibitor Array*

Catalog # IAP002

Instruction Manual

Please read all instructions carefully before use

For Research Use Only

***Patent pending**

www.symansis.com

Introduction

This Easy50™ PI3K Inhibitor Array is a panel of inhibitors that are known to specifically target the PI3K signaling pathway, and related pathways. Data generated with the Easy50™ array produces an inhibition profile that provides a more comprehensive understanding of the specific PI3K signaling pathway involved in a particular response.

These inhibitors are provided in a 96-well plate, and are simply dissolved in DMSO, and then diluted in cell culture medium and added to cells, or diluted in buffer for kinase assays.

Materials provided

1. One 96-well microplate containing inhibitors as shown in the **Plate setup**. The microplate is sealed under argon with aluminium sealing foil.
2. One plate cover to seal plate while inhibitors are dissolving in DMSO.
3. 12 strips of caps to cover the wells of the microplate once inhibitors are dissolved in DMSO.

Materials required but not provided

DMSO (dimethyl sulfoxide), approximately 3.5 ml.

Storage and stability

Store microplate at ambient temperature (<25°C) and use within 1 year. After re-suspension of the inhibitors in DMSO, seal and store at -20° in the foil pouch provided, and use within 3 months. Avoid repeated freeze-thaw cycles. Hygroscopic. Protect from light.

Safety precautions

Avoid contact with skin. Wear appropriate personal protection (labcoat, gloves, safety glasses). For specific safety instructions, refer to the material safety data sheet provided.

Directions for use

see page 5

Data interpretation

To determine the EC₅₀ construct a graph of the cellular response you are measuring versus the diluted drug concentration. The drug concentration at which the response is inhibited 50% is the EC₅₀.

If the EC₅₀ correlates closely with IC₅₀ (see **Inhibitor information** for IC₅₀ values) then this provides evidence that the enzyme targeted by the drug is required for the cellular process you measure.

The diagram below indicates which intermediate proteins are targeted by the various inhibitors. Compare the inhibition profile obtained from a screening with the diagram to identify the probable intermediate proteins in the PI3K pathway involved in a particular response.

Easy50™ PI3K Inhibitor Array




Plate setup



The following table shows the concentrations of each inhibitor after re-suspension in 35 µl of DMSO (= stock solution). The recommended final or working concentration is 1/1000 of the stock solution.

[illegible]



Directions for use



- Open the pouch and remove the microplate. Orientate the plate as shown in **Plate setup** and carefully remove the sealing foil.
- Add 35 μ l of DMSO to each well. (35 μ l will achieve the stock concentrations shown in the **Plate setup**).
- Cover plate with sealing film provided, then gently dissolve inhibitors by placing the plate on a horizontal orbital microplate shaker set at low speed (400 rpm) for 5 minutes. Alternatively, dissolve by pipetting 5x (10-20x for LY294002).
- Remove sealing film, and cap any unused columns using capstrips.
- Pipette appropriate volumes of inhibitor stock solutions into cell culture medium or kinase buffer. A 1/1000 dilution is required (ie. for 2 ml of working solution, add 2 μ l of inhibitor stock solution).
- Ensure that the inhibitor stock solution is fully dissolved in medium / buffer. If the inhibitor solution is pipetted directly into cell culture wells, mix by moving the plate gently to swirl the medium around the well.
- The appropriate incubation time of inhibitors with cells or enzymes should be determined by the investigator depending on their research needs. As a guide, the Easy50 PI3K array inhibitors have been tested with a 10-minute incubation time before stimulus of signaling pathways.
- The bottom well of each column of an inhibitor is empty, and can be used as an experimental control by adding DMSO as for the other wells.
- After use, cap the opened wells with the cap-strips provided, pressing firmly to seal. Return the plate to the foil pouch, and store at -20° .



Inhibitor information



PIK75	CAS Number [372196-67-3]	Chemical Formula: C ₁₆ H ₁₄ BrN ₅ O ₄ S	Molecular Weight: 452.28
	 + + + + + + + + + + 50 100 500 1000 nM Insulin, 1μM, 5mins L6 cells p~Akt^{Ser473} Akt		
Background information PI3K alpha-selective inhibitor with an in vitro IC₅₀ of 8 nM for p110 alpha, 343 nM for p110 delta, and 907 nM for p110 beta¹. 1. Chaussade C, et al. [2007] Biochem J. 404(3):449-58. Evidence for functional redundancy of class IA PI3K isoforms in insulin signalling.			



Note: PIK75, TXG221, PIK294, AS252424, PI103, KU0063794 and 690693 are also available from Symansis in 1-100 mg quantities.



AS-252424	CAS Number [900515-16-4]	Chemical Formula: C ₁₄ H ₈ FNO ₄ S	Molecular Weight: 305.28
			
Background information PI3K gamma-selective inhibitor with an in vitro IC ₅₀ of 30 nM for p110 gamma, 940 nM for p110 alpha, and >20 µM for p110 beta and delta ¹ . 1. Pomel V, et al. (2006) J Med Chem. 49(13):3857-71. Furan-2-ylmethylene thiazolidinediones as novel, potent, and selective inhibitors of phosphoinositide 3-kinase gamma.			

LY-294002	CAS Number [154447-36-6]	Chemical Formula: C ₁₉ H ₁₇ NO ₃	Molecular Weight: 307.34
	 p~Akt^{Ser473} Akt + + + + + + + + + + 0.5 1 5 10 50 µM Insulin, 1µM, 5mins L6 cells		
Background information PI3K inhibitor with an in vitro IC ₅₀ of 500 nM for p110 alpha, 570 nM for p110 delta, and 973 nM for p110 beta ¹ . 1. Chaussade C, et al. (2007) Biochem J. 404(3):449-58. Evidence for functional redundancy of class IA PI3K isoforms in insulin signalling.			

PI-103	CAS Number [371935-74-9]	Chemical Formula: C ₁₉ H ₁₆ N ₄ O ₃	Molecular Weight: 348.36
			
Background information A kinase inhibitor with an in vitro IC ₅₀ of 2 nM for DNA-PK, 8 nM for p110 alpha (PI3K), and 20 nM for mTORC1 ¹ . 1. Knight ZA, et al. (2006) Cell. 125(4):733-47. A pharmacological map of the PI3-K family defines a role for p110alpha in insulin signaling.			

Rapamycin	CAS Number [53123-88-9]	Chemical Formula: C ₅₁ H ₇₉ NO ₁₃	Molecular Weight: 914.17
	 p~p70 S6K^{Thr389} p70 S6K + + + + + Insulin, 1μM, 5mins 50 100 500 1000 nM J774 cells		
Background information Immunosuppressant drug and inhibitor of mammalian-Target-Of-Rapamycin, with an in vitro IC ₅₀ of <10 nM for mTORC1 ¹ . 1. Guertin, DA and Sabatini, DM (2009) Sci. Signal. 2 [76] pe2. The Pharmacology of mTOR inhibition.			

KU0063794	CAS Number [938440-64-3]	Chemical Formula: C ₂₅ H ₃₁ N ₅ O ₄	Molecular Weight: 465.54
	 p~mTOR^{Ser2481} mTOR + + + + + Insulin, 1μM, 5mins 10 50 100 500 nM CHO-IR cells		
Background information Inhibitor of mammalian-Target-Of-Rapamycin complexes, with an in vitro IC ₅₀ of 10 nM for mTORC1 and mTORC2 ¹ . 1. García-Martínez, JM et al. [2009] Biochem. J 421: 29–42. Ku-0063794 is a specific inhibitor of the mammalian target of rapamycin (mTOR).			

690693	CAS Number [937174-76-0]	Chemical Formula: C ₂₁ H ₂₇ N ₇ O ₃	Molecular Weight: 425.48
	 p~GSK3a/b^{Ser21/9} GSK3b + + + + + Insulin, 1μM, 5mins 0.1 0.5 1 5 μM CHO-IR cells		
Background information Pan-Akt kinase inhibitor with an in vitro IC ₅₀ of 2nM for Akt1, 13 nM for Akt2, and 9 nM for Akt3 ¹ . 1. Heering, DA et al. [2008] J. Med. Chem., 51:5663–5679. Identification of 4-[2-[4-Amino-1,2,5-oxadiazol-3-yl]-1-ethyl-7-[[[3S]-3-piperidinylmethyl]oxy]-1H-imidazo[4,5-c]pyridin-4-yl]-2-methyl-3-butyn-2-ol [GSK690693], a novel inhibitor of AKT kinase.			